

Singa® fits you, too.

The Singa® Diaphragm is a contraceptive for all women who want natural and safe contraception.

It is **especially** suitable for women who:

- want or need hormone-free contraception
- cannot or do not want to use an IUD
- are breastfeeding
- want contraception only when needed
- have given birth / are in the interval between two pregnancies
- have a latex allergy
- want combination with natural methods of family planning (e. g. sympto-thermal method)

Singa® can NOT be used:

- within the 6 weeks after childbirth
- in case of marked prolapse of pelvic floor / uterus

Note:

The correct size of the Singa® Diaphragm must be checked by a healthcare professional before purchase. Singa® does not protect against sexually transmitted diseases.

We recommend to always use Singa® in combination with Caya® Gel.

Benefits of Singa® Diaphragm

- hormone-free
- easy to use
- no side effects
- reusable for two years
- latex-free
- ecological
- no impairment of sexual feelings of both partners
- to be applied only when needed
- together with the diaphragm gel it provides **double safety:**
 1. cervical cup as mechanical barrier
 2. cellulose in diaphragm gel as additional mechanical barrier

Together with the Singa® Diaphragm you get:

- instruction for use
- storage case



Singa®

Diaphragm



Caya®

Diaphragm gel



hormone-free
safe
easy to use



Singa®

Diaphragm

Manufacturer:

KESSEL medintim GmbH
Nordendstr. 82 – 84
64546 Mörfelden-Walldorf
Germany

e-mail: export@medintim.de

www.medintim.de

CE 0483
Medical device

Made in Germany



Medintim

Contraception using the Singa® Diaphragm



The Singa® Diaphragm is a new contraceptive barrier device.


The individual size of the Singa® Diaphragm must be checked by a healthcare professional before purchase.



Specifications of the Singa® Diaphragm

Firm insertion edge

- folds compactly for easy insertion
- stable shape during insertion
- supports cervical cup in posterior fornix



Grip dimples (nubs)

- cue for where to hold and squeeze the Diaphragm
- provides enhanced grip when slippery

Cervical cup

- sized to fit a range of cervical positions
- surrounds the cervix




Labelling

- shows product name and size of the Diaphragm

Many women wish for natural and safe contraception – without hormones!

The new Singa® Diaphragm offers reliable and self-determined contraception precisely when you need it. Singa® is inserted into the vagina before intercourse. It covers the cervix completely and is located between the posterior fornix and the niche on the pubic bone. During and after intercourse, Singa® prevents the sperm from entering the uterus. After intercourse, the Singa® Diaphragm must remain in the vagina for at least 6 hours.



Double safety with Singa®

Singa® Diaphragm and Caya® Gel ensure safety:

1. The cervical cup is the mechanical barrier directly in front of the cervix. The solid rim reliably covers the cervix.
2. The Caya® gel contains cellulose, among other things, which forms an additional mechanical barrier directly at the cervix.

Safety, ease of use and wearing comfort of round Diaphragms like Singa® have been demonstrated in clinical studies.

Singa®
Diaphragm

